

Trade in your old shaker and save!

Receive pickup and disposal of your old shaker at no additional cost when you purchase a qualifying Solaris or MaxQ Shaker*

Combine the reliability and technical innovations of our new Thermo Scientific™ Solaris™ Open Air Shakers and the long-lasting stability of our Thermo Scientific™ MaxQ™ Orbital Shakers, designed for demanding daily laboratory use. These versatile shakers are ideal for solubility studies, extraction procedures, protein expression, and many other applications. So for a limited time, you can now save up to \$835 on the purchase of a new shaker when you trade in your old orbital shaker.



NEW Solaris Open Air Shakers

- Accepts loads up to 42 lb. on the 2000 and 84 lb. on the 4000 with single platforms
- Can be used in a variety of standard incubators, warm rooms, environmental chambers and refrigerators
- Durable triple eccentric drive with 10-year warranty
- Programmable features allow for up to 99 unique saved programs for easy set up and reproducibility
- Touchscreen display allows you to clearly view speed and runtimes concurrently, even from a distance
- Historical monitoring of alerts, alarms and interaction logs help provide greater security and reproducibility for your work
- USB and Ethernet connections for monitoring and data output

MaxQ Digital High-Performance (HP) and Standard Shakers

- Designed to help provide the reliable results required for bacteria and yeast cultures
- Durable triple eccentric drive with 10-year warranty
- HP units designed for large or heavy loads to help maximize production
- Weight capacity of up to 71 lbs. for heavy loading with our HP units
- Capacity up to a 6 L flask for high-throughput shaking
- Options available for continuous or timed runs for added flexibility

* Valid while supplies last for orders received between April 1 and December 31, 2020, from residents of the United States (excluding Puerto Rico), and only in connection with the trade-in of orbital shakers. Suppliers of Life Science products are not eligible to participate in this promotion. Determination of eligibility is at the sole discretion of Thermo Fisher. Subject to manufacturer or dealer inspection of the trade-in shaker. Prior to return, customer MUST fulfill certain obligations and complete N206 decontamination form prior to trade-in, and follow instructions for packaging and staging for pickup. Therefore, customer will fully clean the traded-in unit and decontaminate it of all radioactive, biological, toxic and other dangerous or hazardous materials or substances before it can be picked up. Manufacturer will contact customer to arrange pickup time for trade-in shaker within 10 days of receipt of completed form. The traded-in unit MUST be picked up at time of delivering the new unit. Additional charges may apply if trade-in shaker is not available for pickup on pickup date (which is typically when the new unit is also being delivered). Savings based on manufacturer's suggested list price. Subject to dealer participation. Actual savings may vary depending on the type and condition of the trade-in instrument. Manufacturer reserves the right to discontinue or modify this promotion at any time. This offering cannot be combined with other discounts or promotions. Offer void where prohibited, licensed, or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply. Customer acknowledges that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

This promotion is not available to our current or potential healthcare customers. Healthcare customers include any healthcare practitioner, other healthcare provider or any individual or organization authorized to prescribe, dispense, purchase or influence the acquisition or use of medical devices or supplies for clinical use. Healthcare customers include, without limitation, any "Health Care Practitioner" under the Massachusetts Marketing Code of Conduct Law (defined at 105 C.M.R. §907.004) or a "Health Care Provider" under the Vermont Prescribed Products Gift Ban and Disclosure Law (defined at Vt. Stat. Ann. tit. 18, §4631a(8)). By participating in this promotion, you warrant that you are not a current or potential healthcare customer. Thermo Fisher Scientific reserves the right to determine your status in the event of a dispute.

Find out more at fishersci.com/orbital-shaker-upgrade

Contact us today:

In the United States

Order online: fishersci.com

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000



This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for the customers specific use or application. © 2018-2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **BN20180256 0320**